

Nortev About Us

Nortev Ltd., founded in 2008 by two engineers, Declan Moran and Gavan O' Sullivan, is based in Galway, **Ireland**. Nortev is dedicated to developing innovative inhalation devices for animals with the unique design and technology.

Nortev has since grown to become a leader in the design, manufacture, and commercialization of **portable nebulizers** for horses and companion animals, distributing its products worldwide in Europe, North America, Australasia, and the Middle East. Customers include well-known veterinary clinics, stud farms and horse trainers as well as professional and leisure riders.

Nortev is committed to **serving natural needs** of horses and small animals through constant improvement of our Flexineb devices. The respiratory system of each animal is fundamental for their quality of life, which is why we developed the **Flexineb with your animal in mind**.



Nortev Ltd.
GALWAY, Ireland.

Flexineb
by nortev

Flexineb
by nortev

Serving natural needs

**Your horse –
our technology**



NaszKonik.pl



Nortev Ltd.
Unit 14, Westlink Commercial Park,
Oranmore, GALWAY, Ireland, H91 HD62.



Designed and Manufactured by
Nortev Ltd. in Galway, Ireland





The Flexineb E3 is an equine nebulizer. The Flexineb E3 converts liquid solutions which are intended to act in the lungs into a fine mist to be inhaled by the horse.

Features:



Silent



Portable



Flexible & Comfortable



Safe

... easy to use!



Available in 3 sizes to suit your horse's needs.
Choose from multiple colours.

Why Flexineb E3

The horse's lungs make a **fundamental** contribution to the wellbeing and performance of our horses. Unfortunately, they are under constant attack from external influences.

- Environmental factors such as dust, pollen, mould, pollution and weather conditions
- Stressful situations like stabling, travelling, and competition
- Bacteria and viruses in the air
- And many more

Respiratory tract illnesses are the second leading cause of health problems in horses, surpassed only by musculoskeletal conditions. Unlike external injuries to your horse, irritations of their lungs are invisible. We often only react when the first symptoms, such as coughing, appear. **Preventative** measures need to be taken to avoid the development of irreversible damage to the lungs and to safeguard their health and **quality of life**.

Inhalation is an **easy** and non-invasive way to provide your horse with the targeted support they need, whether preventative, supportive or therapeutic. Make Flexineb E3 part of your horse's **daily routine**.

How Flexineb E3 works

1. Fill the Flexineb E3 Cup
2. Turn on the Flexineb E3 Controller
3. The liquid is converted to a **fine mist** (aerosol) by the Flexineb E3 cup through vibration of the mesh/membrane (90 thousand vibrations per second) and released into the Aerosol Holding Chamber.
4. The available aerosol is **actively inhaled** by your horse
5. Fine particles are **delivered** to your horse's lungs
6. Your horse exhales through the exit valve

The cycle repeats with every breath

